

ORGANIC CHEMISTRY

Experiment 1.

Detection of the hydroxyl group in ethanol (iodoform reaction)




Laboratory equipment:

- test tubes in rack,

Chemicals:

- ethanol,
- 2M NaOH,
- Lugol's solution (I₂ in KI)

Pour about 1cm³ of ethanol into a test-tube, add the same volume of NaOH solution and Lugol's solution to obtain bright yellow colour of solution. Heat this mixture at 60 °C, since the bright yellow precipitate of iodoform appears. (The solution has characteristic smell of disinfectant).





Write down the reactions in the table of a report and the names of the products.

Experiment 2.

Determination of aldehyd group (Fehling's reaction)

Fehling reagent (mixture of I and II Fehling solutions) contains Cu²⁺ ions. It reacts with a aldehyde in redox reaction Cu²⁺ ions undergo reduction to Cu⁺ and an aldehyde undergoes oxidation to an acid.



Laboratory equipment:

- test tubes in rack,


Chemicals:

- aldehyde,
- Fehling's solution I (CuSO₄) solution,
- Fehling's solution II (alkaline solution of sodium potassium tartrate),

Mix the solution I and II in 1:1 ratio. Take 1cm³ of this solution and add 1cm³ of an aldehyde solution. Heat up slowly this mixture until boiling. It should appear brick-red precipitate of Cu₂O. Give the description of the experiment and write the equation of the reaction of aldehyde oxidation.

Experiment 3.

Tollens's reaction (Silver Mirror Test)



Laboratory equipment:

- test tubes in rack,
- beaker,

Chemicals:


- aldehyde (formaldehyde),
- 0.05M AgNO₃,
- 2M NH₄OH

Pour 1cm³ AgNO₃ solution into a test-tube and then add as much of NH₄OH solution, as to dissolution of precipitate and then add 1cm³ of aldehyde. Put the test-tube into a water bath. It should appear metallic silver on the walls of the test tube. Write down the reaction.

Experiment 4.

Detection of ketone group (Gunning's reaction)

This reaction is characteristic for methyloketone.



Laboratory equipment:

- test tubes in rack,

Chemicals:

- ketone (acetone),
- 2M NH_4OH
- Lugol's solution.


Add 0.5cm^3 of ketone solution to 1cm^3 of NH_4OH solution and then few droplets of Lugol's solution until yellow colour appears. Heat the mixture until 60°C and leave it for few minutes. It should appear the yellow precipitate of iodoform with characteristic smell.

Write down the reaction's details.

Experiment 5.

Detection of carboxyl group

Detection of COOH group is possible through its transformation to an ester in reaction with alcohol. The esters have characteristic smell of flowers or fruits.



Laboratory equipment:

- test tubes in rack,

Chemicals:

- carboxyl acid,
- ethanol,
- H_2SO_4 concentrated.


Pour 1cm^3 of acid solution into the test-tube, add 1cm^3 of ethanol and few droplets of concentrated sulphuric acid. The mixture should have a characteristic smell of fruits.

Write down the esterification reaction.

Experiment 6.

The saponification reaction

The saponification reaction is the reaction of long chain fatty acids, like Lauric C12, Myristic C14, Decanoic C10, Palmitic C16 and others and sodium hydroxide:



Laboratory equipment:

- 1 test tube,
- water bath,
- gas burner,
- evaporating dish,
- watchglass,
- tripod,
- glass rod.

Chemicals:

- coconut oil,
- 12M NaOH

Put approximately 2 g (1 teaspoon) of coconut oil into the evaporating dish. Pour 2cm^3 of 12M NaOH to the test tube and heat it in the water bath up to approx. 60°C . Place the evaporating dish with fat on the tripod and heat it until melting. Pour warm hydroxide to the oil and mix it with the glass rod until mixture thickens. Pour out the soap on the watchglass and observe.

20..../....	<i>Name, surname:</i>	<i>Assistant signature</i>
WIMiR		
Gr.	<i>Subject:</i> ORGANIC CHEMISTRY	

No	Experiment	Reactions	Observation, conclusions
1.	Detection of hydroxyl group in ethanol (iodoform reaction)		
2.	Detection of aldehyde group (Fehling's reaction)		
3.	Tollens's reaction (Silver Mirror Test)		
4.	Detection of ketone group (Gunning's reaction)		
5.	Detection of carboxyl group		
6.	The saponification reaction		